




**One-step breast reconstruction** with  
**Pectoralis Major intact**

**BRA**XON<sup>®</sup>  
Muscle sparing one-step reconstruction



# the simpler the better

Braxon® is the result of the experience which has developed in the field of tissue bio-engineering joined with that of clinical practice aiming to create a more conservative surgical procedure.

The study of the biomaterial's regenerative capacity applied together with the constantly ambitious challenges in the field of reconstructive surgery has led to the emergence 5 years ago of an innovative surgical technique which signaled a further enhancement in the field of one-step breast reconstruction.



The one-step breast reconstruction, when indicated, has highlighted important benefits for the patient, who recovers her physical integrity in a single operation, together with significant cost savings for the health economy.

Innovations in the field of biomaterial have contributed in a major way to enabling this new operation to substantially improve the cosmetic outcomes of immediate breast reconstruction.

**BRAXON® IS A SYNTHESIS  
BETWEEN THE MOST ADVANCED  
BIOMATERIAL AND THE MOST  
CONSERVATIVE ANATOMICAL  
IMPACT FOR THE PATIENT.**

# Muscle sparing **WHY**

The disinsertion of Pectoralis results in a deficit of muscle function of arm flexion, internal rotation and adduction. Over time the interference in the synergy of the various muscle groups may result in weakness of humeral-scapular articulation with repercussions in movement, pain and also an impaired cosmetic result.

**THE DISINSERTION OF THE MUSCLE CAUSES BLEEDING, POSTOPERATIVE PAIN AND SEROMA.**

The creation of a sub-muscular pocket for the normal procedures, one-step or two-stage breast reconstruction is longer and more complicated and requires postoperative physiotherapy.





BRA<sup>X</sup>ON<sup>®</sup>

BRA<sup>X</sup>ON<sup>®</sup>  
THE MOST PHYSIOLOGICAL  
IMPLANT

Braxon<sup>®</sup> **WHO**  
Ask your patients.

# Above muscular pocket **HOW**

Recent publications have shown a lower ratio of capsular contractions if the breast prosthesis is wrapped in a biological matrix.


Braxon® is a pre-shaped porcine dermis which allows the tailoring of an ADM pocket around the mammary prosthesis and its fixing above the pectoralis muscle, which is kept intact.

Covered by 3 patents, the Braxon® shape perfectly matches the contours of a silicone prosthesis. Having chosen the size and shape of the implant, the customised preparation of the pre-shaped graft to fit around the silicone implant is carried out on a sterile work surface, while the assistant surgeon checks for haemostasis and places the drain in preparation for the insertion of the implant, neatly wrapped in the Braxon® graft, thereby optimizing the operation time.

# Braxon® **WHEN**

In one-step reconstruction after nipple or skin sparing mastectomy. When there is a good vascularised sub-cutaneous layer.

Previous radiotherapy or co-existing medical conditions such as diabetes or connective tissue diseases are contra-indications. The use of monopolar diathermy in the mastectomy skin flap dissection should be reduced to a minimum to prevent skin flap necrosis.



# Technical characteristics

The Braxon® implant is made of porcine dermal extracellular matrix(\*) with a thickness of 0.6 mm.

The exclusive production process has been developed with the aim of generating a completely natural product (not crossed linked) without the presence of any chemical substance which can amplify the inflammatory response and slow the pathway of tissue regeneration. The native proteic structure provides the benefit of immediate bio-availability for incorporation into the host tissue with lower inflammatory responses such as seroma or the red-skin flare phenomenon seen following the implantation of other biomaterials used in breast reconstruction.

## **Natural, artificial and synthetic matrices.**

The **natural** Braxon® matrix is made of native proteic polymers derived from acellularized tissue of porcine dermis (collagen). The organism recognises it as its own, and transforms it into self tissue through the natural regenerative process (remodelling).

**Artificial** meshes derive from natural polymers, but they are chemically modified for reinforcement (e.g.. Cross-linked meshes). They are tolerated by the organism but do not stimulate any regenerative process.

**Synthetic** meshes, made of polymers obtained by chemical synthesis (Polypropylene with metal coatings for example), are chemically and physically tolerated by the organism but they do not stimulate a regenerative process.

**NOWADAYS THE CONCEPT OF BIOCOMPATIBILITY ISN'T ENOUGH TO DEFINE THE EFFECTIVENESS OF A BIOMATERIAL. ITS PERFORMANCE MUST SURPASS PASSIVE TOLERANCE. IT MUST BE ACTIVE, NOT SOLELY INERT, WITH THE AIM OF POWERING THE BIOLOGICAL PROCESS OF GUIDED TISSUE REGENERATION.**

In order to allow Braxon® to be incorporated into the tissues without an amplified inflammatory response, the production process entails the sublimation (freeze drying) in the final phases of the newly created acellular graft, which allows the complete removal of the liquid chemicals used in its preparation through exact use of pressure and temperature. This process makes the product dry, optimal for correct conservation at room temperature, and only requires a simple rehydration before use, without the need for repeat washings to attempt the removal of damaging chemicals which other products require.

Its patented shape allows a perfect fit around the silicone prosthesis which creates a smooth surface with no graft over-lap after simple suturing to the muscle surface. It fits neatly under the skin in a perfect position, checked by the elevation of the patient prior to closing the skin wound.

## **Implantation technique**

Braxon® must be sutured to the pectoralis major with single stitches to ensure primary stability of the matrix which requires intimate contact with vascularized tissue.


## **Postoperative management**

Early and prolonged use of a conforming bra for 3 or 4 weeks as well as a compressive dressing and reduced mobility of the arm will significantly decrease the seroma formation.


## **(\*)Porcine extracellular matrix:**

A complex network of extracellular macromolecules that, in addition to performing a cementing function between cells and tissues, provides an organized structure in which the cells can migrate and interact with each other.

# THE WINNING IDEA



CASE OF BREAST IMPLANT SIZE REVISION:  
IN THIS OCCASION A WELL-VASCULARIZED  
NEO-FASCIA WAS OBSERVED AFTER 7  
MONTHS FROM THE BRAXON<sup>®</sup> PROCEDURE.



# Surgical Steps






During mastectomy surgeon begins Braxon<sup>®</sup> procedure



Braxon<sup>®</sup> must be hydrated for 5 minutes to make it soft and pliable






Use of sizer to check that size and shape of mammary implant is correct



## Tailoring Braxon<sup>®</sup>

Braxon<sup>®</sup> is pre-shaped in such a way as to contain a breast implant of any size and shape. Its use is intuitive and requires scissors and suture to “dress” the prosthesis and be sutured over the pectoralis major muscle.





Surgeon performs Braxon® tailoring around mammary prosthesis



Braxon® is inserted and sutured above the pectoralis major






Skin closure

35'

45'

FINISH  
60'



# Q&A

## Questions & Answers

### Is Braxon® supplied already shaped?

Yes – It is supplied pre-shaped to wrap around different sizes of mammary prosthesis.

### How is the Braxon® customized to fit precisely around the chosen implant?

The surgeon adjusts the pre-shaped Braxon® on a sterile work surface as shown in the brochure.

### Can I use different implant sizes and shapes?

Yes the pre-shaped Braxon® is adjusted to fit around any shape or sized implant.

### What sutures should I use?

Interrupted 3 0 absorbable soft braided sutures.

### How to suture?

Once the tailored customized ADM has been sutured in place around the chosen implant, the superior and upper medial and lateral edges are sutured onto the pectoralis muscle having elevated the patients upper body to check symmetry with the opposite breast.

### How many drains and for how long?

The inflammatory response to the new Braxon® chemical-free ADM is much less, and as a result our experience has shown the drains can be safely removed within a few days, and certainly less than the 2 weeks necessary when alternative products are used.

### What about capsular contraction and cosmesis?

To date, with follow-up at 5 years there has been no clinical evidence of capsular contraction around the sub-cutaneous placement of the Braxon®-wrapped implant.

### Do I need to consider using a round implant to avoid the risk of rotation often seen in sub-muscular implant reconstruction?

No, since the Braxon®-wrapped implant is on top rather than underneath the muscle, it is not subjected to the rotational forces of muscular contraction, and also looks more natural when a shaped implant is employed. Nevertheless Braxon® can perfectly “dress” also a round implant, with same good cosmetic results.



Braxon® bilateral implantation with anatomical implants  
15 months postoperative photographs




# References




ADM Braxon® Pre-shaped for total coverage of the breast implant.


## BRX06S



I.



III.



II.

### Histologies

I. Hematoxylin-eosin staining of the sample implants Braxon® sterile, non-implanted. 10x magnification. The staining reveals the complete absence of cellular material. Observed transverse sections of pre-existing blood vessels which, while retaining their structure, ensure a more rapid permeation of the blood following implantation.

II. Azan-Mallory staining of sample Braxon® sterile prosthesis. 5X magnification. Highlights the collagen fibres of the matrix. Observed absence of cellular material.

III. Hematoxylin-eosin staining of the sample implants Braxon® 4 weeks following implantation. 5X magnification. The hematoxylin colours in violet cellular components showing a high degree of cellular infiltration.

## Bibliography

- Berna G, Cawthorn SJ, Papaccio G, Balestrieri N. Evaluation of a novel breast reconstruction technique using the Braxon® acellular dermal matrix: a new muscle-sparing breast reconstruction. ANZ J Surg, 2014.
- Fahad M. Iqbal, Anjali Bhatnagar, Raghavan Vidya. Host Integration of an Acellular Dermal Matrix: Braxon Mesh in Breast Reconstruction. Clinical Breast Cancer, 2016.
- Maruccia M, Mazzocchi M, Dessy LA, Onesti MG. One-stage breast reconstruction techniques in elderly patients to preserve quality of life. Eur Rev Med Pharmacol Sci, 2016.
- Vidya R, Iqbal FM. A guide to pre-pectoral breast reconstruction: a new dimension to implant based breast reconstruction. Clinical Breast Cancer, 2017.
- Berna G, Cawthorn S. Absence of capsular contracture 4 years after prepectoral breast reconstruction with Braxon® ADM: a case series. European Journal of Plastic Surgery, 2017.
- Vidya R, Masia J, Cawthorn S, et al. Evaluation of the effectiveness of the prepectoral breast reconstruction with Braxon® dermal matrix: first multicentre European report on 100 cases. The Breast Journal, 2017.



Patented by



info@decomed.it



---

[www.braxon.com](http://www.braxon.com)