

## For safe patient transfer and handling

**Air**  
*Slide*



**NEW!**

Reduces manual handling for staff,  
increases comfort for the patient.

The ONLY dual-technology lateral transfer system on the market, combining the Comfort of an air mattress with the quiet, back-saving convenience of a slide sheet. PLUS, RF welding and the absence of air holes, stitching and velcro ensures maximum infection control.

### Key features

- Maximum user weight 544 kg
- No manual lifting
- Protects fragile skin
- Can be used on the floor
- Suitable for use with X-Ray machines
- Extension pull handles to avoid stooping
- Easy to clean
- Portable System - with trolley bag
- No stitching, no air holes, no Velcro, no noise, no weight limit
- Two sizes available - large and bariatric





For safe patient transfer and handling



## Transfer and repositioning are simple and quiet

The AirSlide is an Air Mattress and Slide system all in one device. It can be used with or without air so patient transfer and repositioning are simple and quiet.



### Reduces Injury

The AirSlide system reduces back injuries that occur when transferring a patient from one surface to another. Lateral transfers can be accomplished with fewer personnel while keeping the patient safe and comfortable. With the AirSlide system, medical staff can move a patient with ease throughout the facility for various procedures. With no metal parts, AirSlide even travels safely through X-ray and MRI.

### Infection Control

Infection control reaches a new level. Older technology that utilizes air-flow as a lubricant requires literally thousands of pin holes, thousands of needle and thread stitches and even Velcro. Our patented technology removes all pin holes, stitching and Velcro, so microbes have nowhere to hide! Infections kill over 100'000 people every year. Treating hospital infections costs many Billions a year. M.R.S.A. - and other lesser known infections are able to live in places that are hard to clean, such as pin holes, needle and thread stitching and Velcro.

A washable re-usable cover sheet, or a disposable cover sheet is absorbent while protecting the surface beneath. These keep the mattress clean and in use, rather than desinfect.

### Technical information

|                   | Large 81 cm<br>(32") | Bariatric 117 cm<br>(46") |
|-------------------|----------------------|---------------------------|
| Length            | 198 cm (78")         |                           |
| Max. user weight  | 544 kg (86 stone)    |                           |
| Mattress weight   | 1.6 kg (3.5 lbs)     | 2.2 kg (4.9 lbs)          |
| Mattress height   | 15 cm (6") inflated  |                           |
| Air supply weight | 3.5 kg (7.8 lbs)     |                           |



[www.samarit.com](http://www.samarit.com)

Ask for separate brochures for the other products in the Samarit range!



Samarit Rollbord



Samarit Transglide



Samarit Rollover



Samarit Nursing Aids

Download your language from our website or ask your local distributor.



**Samarit Medizintechnik AG**

Dorfplatz 2 / PO Box 213

CH-8126 Zumikon, Switzerland

Tel.: +41 44 918 10 11 Fax: +41 44 918 24 39

e-mail: [info@samarit.com](mailto:info@samarit.com) - [www.samarit.com](http://www.samarit.com)