



□ www.edafico.ro

Sibedafico Medical srl

Str. D.D. Rosca nr.70,
Sibiu, Jud.Sibiu
Fax: + 40 378.106.207
CUI: RO 32732861
Nr. ONRC J32/103/2014
BRD Sibiu
RO54BRDE330SV75516613300

1. IMBRACAMINTE PENTRU BLOCUL OPERATOR

Halat chirurgical (SMS), steril

Captusit :

Cod articol	mar.	buc./cutie	pret unitar
6621481	L	40	15.00 RON
6621581	XL	40	15.00 RON
6621681	XXL	40	15.00 RON



Necaptusit :

Cod articol	mar.	Buc./cutie	pret unitar
6621480	L	40	14.00 RON
6621580	XL	40	14.00 RON
6621680	XXL	40	14.00 RON



2. SETURI SI CAMPURI DE ACOPERIRE PENTRU BLOCUL OPERATOR

2.1 CHIRURGIE GENERALA

SET UNIVERSAL


PRET UNITAR : 72.00 RON

Cod articol : 56 01 001

Cant.	Denumire	Marime in cm
1	Camp chirurgical cu banda adeziva, tristratificat	150 x 230
1	Camp chirurgical cu banda adeziva, tristratificat	175 x 180
2	Campuri chirurgicale cu banda adeziva, tristratificat	75 x 100
1	Banda adeziva	9 x 50
1	Sac protectie masa instrumente	80 x 140
4	Servetele	20 x 25
1	Camp chirurgical bistratificat	75 x 90

Impachetat in :

1	Camp chirurgical pentru impachetare	140 x 190
1	Punga sterilizare	
Nr. buc. / cutie : 12		



SET UNIVERSAL PEDIATRIC


PRET UNITAR: 58.00 RON

Cod articol : 56 01 008

Cant.	Denumire	Marime in cm
1	Camp chirurgical cu banda adeziva, tristratificat	150 x 160
1	Camp chirurgical cu banda adeziva, tristratificat	180 x 175
2	Campuri chirurgicale cu banda adeziva, bistratificat	50 x 50
1	Camp chirurgical, bistratificat	75 x 90
1	Banda adeziva	9 x 50
1	Sac protectie masa instrumente	80 x 140
4	Servetele	20 x 25

Impachetat in :

1	Camp chirurgical pentru Impachetare	140 x 190
1	punga sterilizare	
Nr. buc. / cutie : 10		



2.2 . ORTOPEDIE / TRAUMATOLOGIE

SET ARTROSCOPIE


PRET UNITAR: 95.00 RON

Cod articol : 56 02 004

Cant.	Denumire	Marime in cm
1	Camp chirurgical extremitati, tristratificat, Orificiu elastic central : 5 x 8 cm	235 x 320
2	Benzi adezive	9 x 50
1	Saculet ortopedic	38 x 72
1	Sac protectie masa intrumente	80 x 140

Impachetat in :

1	Hartie sterilizare
1	Punga sterilizare
Nr. buc. / cutie : 14	



SET ARTROSCOPIE „HIGH PROFESSIONAL”


PRET UNITAR: 118.00 RON

Cod articol : 56 02 005

Cant.	Denumire	Marime in cm
1	Camp chirurgical extremitati, tristratificat, Orificiu elastic central : 5 x 8 cm	235 x 320
1	Camp chirurgical, bistratificat	90 x 150
1	Saculet ortopedic	38 x 72
1	Banda adeziva	9 x 50
1	Scai 25	
1	Sac protectie aparat, cu ajutor de introducere	13 x 250
1	Sac protectie masa instrumente	80 x 140
1	Halat chirurgical, SMS, captusit XL	

Impachetat in :

1 Hartie sterilizare
1 Punga sterilizare
Nr. buc. / cutie : 10



SET ARTROSCOPIE CU PUNGA COLECTOARE


PRET UNITAR: 128.00 RON

Cod articol: 56 02 013

Cant.	Denumire	Marime in cm
1	Camp chirurgical ptr.artroscopie orificiu central elastic 5 x 8 cm punga colectoare	235 x 310
1	Camp chirurgical bistratificat	90 x 150
1	Sac ortopedic bistratificat	37.5 x 72
1	Banda adeziva elastica	9 x 50
1	Scai 25	
1	Husa protectie aparat	13 x 250
1	Sac masa instrumente	80 x 140

Impachetat in:

1 Camp de ambalare 140 x 190
1 Punga sterilizare
Nr. buc / cutie : 8



SET SOLD

PRET UNITAR: 148.00 RON

Cod articol : 56 02 001


Cant.	Denumire	Marime in cm
1	Camp chirurgical cu slit si banda adeziva, tristratificat, Slit : 6 x 100 cm	200 x 250
1	Camp chirurgical cu banda adeziva, tristratificat	180 x 290
1	Camp chirurgical, bistratificat	150 x 180
1	Camp chirurgical, bistratificat	75 x 90
2	Benzi adezive	9 x 50
1	Saculet ortopedic, mare	30 x 100
1	Sac protectie masa instrumente	80 x 140
4	Servetele	20 x 25

Impachetat in :

1 Camp de ambalare 140 x 190

1 Punga sterilizare

Nr. buc / cutie : 6



SET EXTREMITATI

PRET UNITAR: 125.00 RON

Cod articol: 56 02 010


Cant.	Denumire	Marime in cm
1	Camp chirurgical ptr.artroscopie orificiu central elastic 5 x 8 cm punga colectoare	235 x 310
1	Camp chirurgical bistratificat	90 x 150
1	Camp chirurgical adeziv, bistratificat	75 x 75
2	Halat standard XL, SMS	
2	Banda adeziva elastica	9 x 50
1	Sac masa instrumente	80 x 140
4	Servetele	20 x 25

Impachetat in:

1 Camp de ambalare 140 x 190

1 Punga sterilizare

Nr. buc / cutie : 8



2.3 OFTALMOLOGIE

SET OFTALMOLOGIE

PRET UNITAR: 55.00 RON

Cod articol : 56 03 200

Cant. Denumire

- | | |
|---|---|
| 1 | Camp chirurgical oftalmologic din material Spunlace cu folie de incizie Integrata, orificiu : 7 x 10 cm si 1 puna oftalmologica |
| 4 | Servetele |

Marime in cm

120 x 150


20 x 25 cm

Impachetat in :

1 Camp de ambalare 110 x 140

1 Punga sterilizare

Nr. buc./ cutie : 40



CAMP OFTALOMOLOGIC

PRET: 24.50 RON

Cod articol : 55 03 202

Cant. Denumire

- | | |
|---|--|
| 1 | Camp chirurgical oftalmologic bistratificat orificiu adeziv si puna colectoare |
|---|--|

Marime in cm


150 x 105

Impachetat in :

1 Hartie sterilizare

1 Punga sterilizare

Nr. buc./ cutie : 50



2.4. CHIRURGIE CARDIACA

Pentru chirurgia cardiovasculara va oferim sisteme speciale si complete de izolare. Dupa alegere, numai un singur camp de acoperire sau intreg sistemul.

SET ANGIOGRAFIE

PRET UNITAR : 99.00 RON

Cod articol : 56 05 502

Continut :

1 camp angiografic, tristratificat, 120 x 300 cm,
cu doua orificii adezive de 9 cm, cu 2 folii de 70cm pe ambele parti.

1 husa tub, 120 x 60 cm, deschidere elastica

1 sac polietilena, 110 x 80 cm, deschidere elastica

1 husa consola, 90 x 60 cm, deschidere elastica

1 halat chirurgical, marimea XL, captusit,


2 servetele 20 x 25 cm

ambalat in :

1 camp chirurgical, 140 x 110 cm

1 punga pentru sterilizare

Nr.buc/cutie: 10



SET CHIRURGICAL ANGIOGRAFIC

PRET UNITAR: 118.00 RON


Cod produs:5605504

Continut:

- 1 camp angiografic, albastru tristratificat, 120 x 300 cm, cu doua orificii adezive de 9 cm, cu 2 folii de 70 cm pe ambele parti.
- 1 husa consola, 120 x 120 cm, cu deschidere elastica
- 1 sac polietilena, 85 x 90 cm, cu deschidere elastica
- 1 husa consola, 90 x 60 cm, cu deschidere elastica
- 1 halat chirurgical, marimea XL, captusit,
- 2 halat chirurgical, marimea XL, standard
- 1 servetel 20 x 25 cm

Impachetat in :

- 1 camp chirurgical, 140 x 110 cm
 - 1 punga pentru sterilizare
- Nr.buc/cutie: 8



SET CARDIOLOGIE PEDIATRICA

PRET UNITAR: 205.00 RON

Cod produs : 55 05 104

Componenta :


- 1 Camp chirurgical tristratificat, 235 x 275 cm, orificiu decentral adeziv 8 x 15 cm
- 1 Camp chirurgical adeziv, tristratificat, 100 x 100 cm
- 4 Camp chirurgical adeziv, tristratificat, 180 x 150 cm
- 3 Banda adeziva elastica 9 x 50 cm
- 3 Punga colectoare autoadeziva 33 x 24 cm
- 1 Halat chirurgical captusit mar. L, SMS
- 2 Halat chirurgical captusit mar. XL, SMS
- 1 Halat chirurgical captusit mar. XXL, SMS

Impachetat in :

1 camp acoperire masa 140 x 190 cm

1 punga sterilizare

Nr.buc/cutie : 4



SIBEDAFICO MEDICAL


SET CARDIOLOGIE PEDIATRICA

PRET UNITAR : 129.00 RON

Cod produs : 55 05 103

Componenta :

1 Camp chirurgical tristratificat, 235 x 275 cm,
orificiu decentral adeziv 6 X 10 cm
1 Camp chirurgical adeziv, tristratificat, 100 x 100 cm
3 Camp chirurgical adeziv, tristratificat, 180 x 150 cm
1 Camp chirurgical tristratificat, 150 x 180 cm
3 Banda adeziva elastica 9 x 50 cm
3 Punga colectoare autoadeziva 33 x 24 cm
1 Halat chirurgical captusit mar. L, SMS
2 Halat chirurgical captusit mar. XL, SMS
1 Halat chirurgical captusit mar. XXL, SMS
Impachetat in :
1 camp acoperire masa 140 x 190 cm
1 punga sterilizare
Nr.buc/cutie: 4



SET PACE MAKER

PRET UNITAR: 122.00 RON

Cod produs : 56 05 404

Componenta :

1 Camp chirurgical tristratificat 90 x 360 cm, orificiu adeziv decentral 10 x 13 cm, folie transparenta de 70 cm

1 Husa folie transparenta 100 x 90 cm, inchidere elastica

1 Caserola rotunda 500 ml

2 caserole rotunde 250 ml

2 betigas betadina

2 servetele 20 x 25 cm

2 Banda adeziva elastica 9 x 50 cm

1 Camp acoperire masa 140 x 190 cm


1 Halat chirurgical captusit mar.L, Spunlace – ambalat individual in hartie sterilizare

Impachetat in :

1 Camp acoperire masa 140 x 190 cm

1 Punga sterilizare

Nr.buc/cutie : 4



SET ANGIOGRAFIE CU COMPONENTE

PRET UNITAR: 116.00 RON

Cod produs : 56 05 511

Componenta :


- 1 Camp chirurgical angiografic 150 x 290 cm, SMS,
- 2 orificii femurale adezive de 12 cm, cu folie transparenta de 70 cm
- 1 Camp chirurgical bistratificat 50 x 50 cm
- 1 Boneta folie transparenta 120 x 120 cm, inchidere elastica
- 1 caserola antireflux
- 1 Caserola rotunda 500 ml
- 2 caserole rotunde 250 ml
- 1 tava 41.3 x 37.2 x 2.7 cm

Impachetat in :

- 1 Camp acoperire masa 140 x 110 cm

- 1 Punga sterilizare

Nr.buc/cutie : 3



2.5.GINECOLOGIE

In special pentru aceasta specializare va oferim o gama completat de sisteme de izolare. Va puteti alege din modele noastre conform metodei dvs. de operare.

SET CHIRURGICAL CEZARIANA

PRET UNITAR: 131.00 RON


Cod articol : 56 06 100

Cant.	Denumire	Marime in cm
1	Camp chirurgical pt. Cezariana, tristratificat, Orificiu adeziv , saculet integrat	235 x 320
2	Banda adeziva elastica	9 x 50
1	Sac protectie masa instrumente	80 x 140
4	Servetele	20 x 25

Impachetat in :

1	Camp de ambalare	140 x 190
1	Punga sterilizare	

Nr. buc./ cutie : 10



SET CHIRURGICAL CEZARIANA

PRET UNITAR : 128.00 RON


Cod articol : 56 06 101

Cant.	Denumire	Marime in cm
1	Camp chirurgical pt. Cezariana, tristratificat, Orificiu adeziv , saculet integrat	235 x 310
1	Banda adeziva elastica	9 x 50
1	Sac protectie masa instrumente	80 x 140

Impachetat in :

1	Camp de ambalare	140 x 190
1	Punga sterilizare	

Nr. buc./ cutie : 10



SET GINECOLOGIE – CISTOSCOPIE

PRET UNITAR: 44.00 RON


Cod articol : 56 06 206

Cant.	Denumire	Marime in cm
1	Camp chirurgical ginecologic, tristratificat, Orificiu : 16 x 20 cm	75 x 200 cm
2	Saculete	50 x 110 cm
1	Camp acoperire	75 x 90 cm

Impachetat in :

1	Camp de ambalare	140 x 190 cm
1	Punga sterilizare	

Nr. buc./ cutie :30



SET ASISTENTA NASTERE

PRET UNITAR: 71.00 RON

Cod articol : 56 06 207

Continut:

6 Campuri chirurgicale autoadezive, tristratificate

1 sac protectie masa instrumente

1 halat chirurgical XL, captusit SMS

Marime in cm

75 x 90 cm


80 x 140 cm

impachetat in :

1 hartie ptr sterilizare

1 punga sterilizare

Nr. buc. / cutie: 10



2.6. UROLOGIE

SET UROLOGIE

PRET UNITAR : 112.00 RON


Cod articol : 56 06 300

Cant.	Denumire	Marime in cm
1	Camp chirurgical urologic, tristratificat, Orificiu : 5 / 9 cm, 20 / 20 cm cu plasa, degetel, sita si saculet integrat	235 x 250 cm
1	Camp chirurgical de acoperire, bistratificat	90 x 150 cm
4	Servetele	20 x 25 cm

Impachetat in :

1	Camp de ambalare	140 x 190 cm
1	Punga sterilizare	

Nr.buc./ cutie : 12



SET UROLOGIE

PRET UNITAR : 54.00 RON


Cod articol : 56 06 301

Cant.	Denumire	Marime in cm
1	Camp chirurgical urologic, tristratificat, Orificiu adeziv central 8 x 8 cm	75 x 200
2	Saculet ortopedic	50 x 110
1	Banda adeziva	9 x 50
1	Camp chirurgical bistratificat	75 x 90

Impachetat in :

1	Hartie sterilizare
1	Punga sterilizare

Nr.buc./ cutie : 30





3. CAMPURI DE ACOPERIRE SUPLIMENTARE

3.1. CAMPURI CHIRURGICALE

Cod articol	Observatii	Dimensiuni in cm	Nr.buc /cutie	Pret unitar/ buc
550 45 75	bistratificat	45 x 75	400	3.75 RON
550 50 50	bistratificat	50 x 50	400	3.90 RON
550 75 92	bistratificat	75 x 90	200	4.25 RON
550 90 15	bistratificat	90 x 150	120	12.90 RON
550 15 18	bistratificat	150 x 180	50	15.50 RON

3.2. CAMPURI CHIRURGICALE ADEZIVE

Cod articol	Observatii	Dimensiuni in cm	Nr.buc /cutie	Pret unitar / buc
553 45 75	bistratificat	45 x 75	120	4.30 RON
553 50 50	bistratificat	50 x 50	200	4.55 RON
553 75 90	bistratificat 7	5 x 90	160	5.90 RON



3.3. CAMPURI CHIRURGICALE CU ORIFICIU SI BANDA ADEZIVA

Cod articol	Observatii	Dimensiuni in cm	Nr.buc /cutie	Pret unitar / buc
552 75 90	bistratificat, orificiu : Ø 6 cm	75 x 90	100	7.05 RON
552 50 60	bistratificat, orificiu : Ø 6 cm	50 x 60	240	6.05 RON
552 90 90	tristratificat, orificiu : Ø 10 cm	90 x 90	100	8.00 RON

4. ACCESORII

4.1. DIVERSE

Sac masa instrumente

Cod articol	Observatii	Dimensiuni in cm	Nr.buc /cutie	Pret unitar / buc
660 77 14		80 x 140	100	18.00 RON



Saculet ortopedic cu 2 benzi adezive

Cod articol	Observatii	Dimensiuni in cm	Nr.buc /cutie	Pret unitar / buc
550 20 07	bistratificat	38 x 72	60	18.00 RON

4.2. PRODUSE PENTRU PROTECTIA APARATELOR

Husa protectie aparat cu ajutor de introducere

Cod articol	Observatii	Dimensiuni in cm	Nr.buc /cutie	Pret unitar / buc
660 13 24	Pliu telescopic, Varf ascutit, 1 banda adeziva	13 x 243	200	20.00 RON
660 15 25	Pliu telescopic, varf bont, 1 banda adeziva	13 x 250	200	20.00 RON



PRODUSE NESTERILE

HUSA TARGA 100 X 210 CM

PRET UNITAR : 11.00 RON

Cod produs : 332 10 21
Spunbond bistratificat
Nr.buc/cutie:80

COMBINEZON UNICA FOLOSINTA

PRET UNITAR : 7.50 RON

Cod produs: 441 20 00
Nr.buc/cutie : 50

PIJAMA UNICA FOLOSINTA

PRET UNITAR : 14.50 RON

Cod produs : 441 40 06
Nr.buc/cutie : 50

HALAT VIZITATOR

PRET UNITAR : 3.3 RON

Cod produs: 441 62 05
Nr.buc/cutie : 100

**ORICE COMANDA TREBUIE SA RESPECTE MODUL DE AMBALARE PE
CUTIE SAU CARTON**