

News letter



- Sensibilitate 93%

- Specificitate 98%

- Precizie 96% (comparativ cu coloratia Gram/ Scorul Nugent)

- Valoare predictiva pozitiva 90.5%

- Valoare predictiva negativa 97.6%



1994 – 2014 | profesionalism
experiență
inovatie



BVBLUE® Test Rapid pentru Vaginite Bacteriene

- Test rapid de detectie a vaginitelor produse de
 - ✓ *Gardnerella vaginalis*
 - ✓ *Bacteroides spp.*
 - ✓ *Prevotella spp.*
 - ✓ *Mobiluncus spp.*
- Test rapid ce detecteaza activitatea **sialidazei**, produsa de patogeni bacterieni vaginali, in probele de lichid vaginal
- Test rapid foarte precis atat la pacientele simptomatice cat si asimptomatice
- Pacientele pot fi testate rapid in cabinetul dumneavoastra

**Stiati ca activitatea crescuta a
sialidazei masurata la 12 saptamani
de sarcina este puternic asociata cu
nasteri premature spontane precoce
si avorturi tardive?***

*Cauci S, Culhane JF. High sialidase levels increase preterm birth risk among women who are bacterial vaginosis-positive in early gestation. *American J Obstet Gynecol* 2011; 204:142.e 1-9.

Oferta speciala: 25% discount pentru 50 de teste rapide BVBlue*



- Colectare proba lichid vaginal prin metoda standard
- Timp de lucru < 1 minut
- Rezultat in 10 minute, la 17-37°C
- Schimbare de culoare evidenta



Denumire	Cod produs	Configurare	Pret RON fara TVA
Teste Rapide BVBlue	BVB25	50 teste + Control kit	400**
Kit control BVBlue	S751		

*in limita stocului disponibil, pana la 30.09.2014

**8 ron fara TVA /testare

SC Novaintermed SRL

Strada Drumul Potcoavei nr. 5A, Pipera-Tunari 077190, Judetul Ilfov
Tel.: +4 0740 309 406 ; Fax: +4 031 401 10 89

sanatateafemeii@novaintermed.ro

www.novaintermed.ro