

News letter



**14,3
Ron/test
(faraTVA)**



Test imunocromatografic in-vitro, calitativ, vizual si rapid, pentru detectia diferentiata a anticorpilor IgG si IgM la virusul Zika.

Test rapid IgG/IgM pentru virusul Zika (A03-34-322)

Metoda de detectie- Imunocromatografica	1 kit =25 teste
Tip proba: ser/plasma si sange integral	Valabilitate kit-minim 18 luni
Conditii de stocare: 2-30°C	Timp citire rezultat: 10 min

Continut pachet: caseta de test, desicant
Tampon pentru proba (3 ml), flacon pentru 25 teste.

Dispozitivul de testare contine un control intern de calitate.
*Se recomanda utilizarea zilnica a unui control extern pentru a fi asigurate performantele optime ale dispozitivului.

Marcaj :CE

