



OZONETTE

medical ozone generator

The OZONETTE is a compact medical ozone generator intended for those applications where safe and accurate concentrations are needed. Measures are shown in standard normalized units. OZONETTE has been designed bearing in mind needs of those healthcare professionals focused on some specific medical applications as described.

Further to this, OZONETTE provides the manual versatile operation mode for applications such as Major Autohemotherapy (MAHT), allowing to set ozone dosage and concentration independently one from the other. Infiltrations are made simple by automatic syringe detection without needing previous flushing because OZONETTE has no gas dead volume, allowing faster and safer fillings of wide range of syringe sizes.

OZONETTE has a self monitoring system showing an error message when the unit is out of specifications providing the practitioner safer and reliable applications.



MEDICAL APPLICATIONS

- Pain Management
- Intra-Articular Infiltrations
- Insufflations
- Bagging Therapy
- Wound Healing



SPECIFICATIONS

Concentration Range	1 to 80ug/Nml in 1ug/Nml steps
Concentration measurement	ACM Maintenance free system
Operation modes	Syringe, Continuous and Manual
Performance	3g/h specially useful for bagging application
Ozone eliminator	Two internal non organic catalyser
Ozone output	Stainless steel Luer lock
Ozone output flow	Digitally adjustable from 10 to 50L/h
Ozone output system	Two valve digitally controlled
Control and operation	Switch less touch pad
Display	High contrast graphic LCD, 192mm x 64mm
Weight / Dimensions	3kg / 25 x 27 x 12cm
Certificates	CE MDD

Ozt_V3_12_15.ing

Super Compact Size:



SEDECAL
 C/ Pelaya, 9 - Pol. Ind. "Río de Janeiro"
 28110 - Algete (Madrid-SPAIN)
 Central Offices:
 Tel: (34) 91 - 628 0544
 Fax: (34) 902 190 385
 Email: O3@sedecal.com

